iba-worldwide.com

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
THIS DOCUMENT IS THE PROPERTY OF I.B.A AND SHALL NOT BE ALTERED COPIED OR COMMUNICATED WITHOUT THE WRITTEN CONSENT OF I.B.A.

S.A. ION BEAM APPLICATIONS

REV:	DATE:	MODIFICATION:	DRAFTMAN:	CHECKED BY:	VALIDATED BY:
0	29/02/16	Original issue	ADUVIVI	QBA	PV
A	01/09/17	New Major Revision	BSOL	QBA	PV
B	26/10/18	New Major Revision	BSOL	QBA	VVA
B.1	07/10/19	See Fig. 00.00-1	ADUVIVI	QBA	VVA

MATERIAL: --
SCALE: --
DIMENSIONS: --
TOLERANCES: --

PROJECT: PROTON THERAPY

S/PROJECT: PROTEUS ONE STD (METRIC)


I. GENERAL INTRODUCTION


TITLE:

Facility General Display 3D Views

07.47.02.

00.00-2 B.1 2/3

REV:	DATE:	MODIFICATION:	DRAFTMAN:	CHECKED BY:	VALIDATED BY:
0	29/02/16	Original issue	ADUVIVI	QBA	PV
A	01/09/17	New Major Revision	BSOL	QBA	PV
B	26/10/18	New Major Revision	BSOL	QBA	VVA
B.1	07/10/19	See Fig. 00.00-1	ADUVIVI	QBA	VVA

THIS DOCUMENT IS THE PROPERTY OF I.B.A AND SHALL NOT BE ALTERED COPIED OR COMMUNICATED WITHOUT THE WRITTEN CONSENT OF I.B.A.

S.A. ION BEAM APPLICATIONS

MATERIAL: --
SCALE: --
DIMENSIONS: --
TOLERANCES: --

PROJECT: PROTON THERAPY

S/PROJECT: PROTEUS ONE STD (METRIC)

I.
GENERAL INTRODUCTION

TITLE:

Facility General Display
3D Views

07.47.02.

00.00-3 B.1 3/3

As a general orientation reference, this figures shows the treatment level floor plan of the facility, with the proton therapy equipment installed. As a convention, we consider the North as indicated on this figure. Everything in this document (text and figures) is coordinated to this reference.


1. The Proton Therapy Equipment (PTE) is installed on three levels: The middle level (Treatment Level) contains most of the PTE (such as the synchro-cyclotron, the beam line, the Treatment Room, Treatment Control Room, PTEV offices). See other levels on next figures.


2. As a building standard, D/BT shall follow normal concrete construction tolerance of +/- 1cm (1/2") unless specifically mentioned in other parts of the text or figures.

3. Caution :It is mandatory to have the Treatment Level 1250 mm under Median Plan level, within a tolerance of +10 -0 mm.
In this IBD, the Treatment Level consists of concrete finished floor. For other ground cover type or floor finish, D/BT shall adapt the floor levels to ensure that the finished floor level complies with this dimension.

4. Shielding thickness is under D/BT responsibility (Indicative value). See ***Fig. 62.01*** for "Sources of Radiation and Shielding".

5. Other PTEV Rooms (not shown on the drawing) shall be located at this level, see ***Figures 35.01 and 36.01***.

TREATMENT LEVEL - FLOOR

- Raised floor
- Grating

TREATMENT LEVEL - CEILING

- Height under concrete 3500 mm
- Height under concrete 5150 mm
- Height under concrete 3800 mm

TREATMENT LEVEL - SUSPENDED CEILING

- Height under ceiling 2500 mm
- Height under ceiling 3500 mm

0 5 10 METER

THIS DOCUMENT IS THE PROPERTY OF I.B.A AND SHALL NOT BE ALTERED COPIED OR COMMUNICATED WITHOUT THE WRITTEN CONSENT OF I.B.A.

S.A. ION BEAM APPLICATIONS

REV:	DATE:	MODIFICATION:	DRAFTMAN:	CHECKED BY:	VALIDATED BY:
0	29/02/16	Original issue	ADUVIVI	QBA	PV
A	01/09/17	New Major Revision	BSOL	QBA	PV
B	26/10/18	New Major Revision	BSOL	QBA	VVA
B.1	07/10/19	See Fig. 00.31-2&3	ADUVIVI	QBA	VVA

MATERIAL: --
SCALE: 1/200 (A3)
DIMENSIONS: mm
TOLERANCES: --

PROJECT: PROTON THERAPY

S/PROJECT: PROTEUS ONE STD (METRIC)

I. GENERAL INTRODUCTION

TITLE:


Facility General Layout Treatment Level

07.47.02.


00.31-1 B.1 1/3

1. The highest level (upper level) includes technical rooms such as the Power Supply Room, the Water Cooling Room, Server Room, Storage Rooms,... ***See Chapter II Other Rooms***.

2. The lower level consists of the pit where the compact gantry structures are supported.

UPPER LEVEL

BASEMENT

Raised floor

Grating (trenches)

Cyclo and CG Hatches

0 5 10 METER

THIS DOCUMENT IS THE PROPERTY OF I.B.A AND SHALL NOT BE ALTERED COPIED OR COMMUNICATED WITHOUT THE WRITTEN CONSENT OF I.B.A.

S.A. ION BEAM APPLICATIONS

REV:	DATE:	MODIFICATION:	DRAFTMAN:	CHECKED BY:	VALIDATED BY:
0	29/02/16	Original issue	ADUVIVI	QBA	PV
A	01/09/17	New Major Revision	BSOL	QBA	PV
B	26/10/18	New Major Revision	BSOL	QBA	VVA
B.1	07/10/19	Upper Level optimization	ADUVIVI	QBA	VVA

MATERIAL: --
SCALE: 1/200 (A3)
DIMENSIONS: mm
TOLERANCES: --

PROJECT: PROTON THERAPY

S/PROJECT: PROTEUS ONE STD (METRIC)

I.
GENERAL INTRODUCTION


TITLE:

Facility General Layout
Upper Level
& Basement


07.47.02.

N


00.31-2 B.1 2/3


SECTION A-A

SECTION C-C

SECTION B-B

SECTION D-D

0 5 10 METER

THIS DOCUMENT IS THE PROPERTY OF I.B.A AND SHALL NOT BE ALTERED COPIED OR COMMUNICATED WITHOUT THE WRITTEN CONSENT OF I.B.A.

S.A. ION BEAM APPLICATIONS

REV:	DATE:	MODIFICATION:	DRAFTMAN:	CHECKED BY:	VALIDATED BY:
0	29/02/16	Original issue	ADUVIVI	QBA	PV
A	01/09/17	New Major Revision	BSOL	QBA	PV
B	26/10/18	New Major Revision	BSOL	QBA	VVA
B.1	07/10/19	Upper Level optimization	ADUVIVI	QBA	VVA

MATERIAL: --
SCALE: 1/200 (A3)
DIMENSIONS: mm
TOLERANCES: --

PROJECT: PROTON THERAPY

S/PROJECT: PROTEUS ONE STD (METRIC)

I.
GENERAL INTRODUCTION

TITLE:

Facility General Layout Sections

07.47.02.

00.31-3 B.1 3/3

Contact us

IBA - Building Integration
Vincent Vaerman
Phone: +32 10 203 427
Email: vincent.vaerman@iba-group.com

www.iba-proteusone.com

iba-worldwide.com

IBA Group | IBAL01/6037 | Proteus®ONE | July 2018

IBA: The best in proton therapy today and tomorrow

Together with our clinical partners, we brought proton therapy to clinical cancer care.

Ever since we started more than 30 years ago, our collaborations, our visionary roadmap and progressively unrivalled experience have enabled us to continue to innovate. Care givers now benefit from leading proton therapy technologies.

Today, our true continuum of Image-Guided Intensity Modulated Proton Therapy solutions can easily be integrated in most healthcare settings to make it available to all patients who need it.

Backed by IBA's unique service offer (financing, workflow optimization, education), our single-room PROTEUS®**ONE**, all our solutions and robust processes (installation, operations and upgrades) are developed in collaboration with our end-users.

Tomorrow, our unique and open culture of sharing will further strengthen the clinical and patient communities we have always cared for. Working collectively, we will achieve our goal which is to offer cancer patients access to effective treatments with decreased side effects and better quality of life.

Request more information: pone@iba-group.com

Visit us online at:

www.iba-protontherapy.com

PROTECT +
ENHANCE +
SAVE LIVES